



forest dental

affordable excellence + designer friendly

Doctor's Cart 5200

\$5,950.00

STANDARD FEATURES

- Standard Control Head
 - automatic control for three handpieces
 - asepsis 3-way air/water syringe with tubing
 - asepsis 4-hole Midwest handpiece tubing
 - pneumatic master control block with diaphragm
 - individual handpiece pressure adjustment
 - non-retraction water coolant with on/off
 - individual water coolant adjustment
 - air coolant adjustment
 - master on/off
 - handpiece flush
 - handpiece exhaust filter/oil collector
 - self-contained clean water system with 1.5 litre bottle
 - wet/dry foot control
- Standard utility centre with housing and cover
- Adjustable height frame from 66cm to 88cm
- Twin wheel hard floor castors
- External 1.1 metre umbilical
- Air pressure regulator with filter
- Master air shut-off valve



AFFORDABLE EXCELLENCE

Forest Dental offers a range of equipment to suit every budget

To talk to one of Forest's specialists call Chris Allen on **0405 813 027**

PO Box 126 Blackheath, Beautiful Blue Mountains, NSW 2785 Australia

Phone: +612 4787 1455 • Fax: +612 4787 1117

Email: info@forestdental.com.au • Web: www.forestdental.com.au